

HotelNewsNow.com

Copyright © 2004-2011 Smith Travel Research /DBA HotelNewsNow.com (HNN). All Rights Reserved. 16 November 2011

Bedbug monitor



The BB Alert Passive Bed Bug Monitor from USBedBugs.com features an optimized bedbug harbourage that encourages them to interact with the monitor. An eco-friendly device, it is heat, chemical, fragrance and upkeep free. Design allows easy inspection for fecal traces, cast skins and live samples in a few seconds. The monitor can be placed on the head end of the bed box springs (if available) between the headboard supports, on the underside of the top slat of the bed at the head end, between the mattress and the bed box springs or platform at the headboard end of the bed, on the rear of the headboard or on the wall below the headboard of the bed, level with the mattress. The monitor should be inspected at least once a month for general monitoring and as often as once a week in high-risk situations. The monitors cover includes a calendar grid that should be used to note the inspection date.